



Pastel Tints

FEATURES




- tints available on:
38%
GMA
60%
67%
SiH
- opaque tints available on:
67%
- clear / open pupils
- black / closed pupils
- standard pupil size 4.5mm




BENEFITS






- colour enhancement of light coloured eyes
- can be applied to both iris and pupil
- opaque backing on iris will create a cosmetic lens which will change the iris colour









Reactive dye tints available on all Hydrocyl, Hydrofit and Hydrolens designs in selected materials.







Tint Range




Apple Green 1	Apple Green 2	Apple Green 3
		
38, 60, 67, SiH	38, 60, 67, SiH	38, 67, SiH

Aqua 1	Aqua 2	Aqua 3
		
38, 60, 67, SiH	38, 60, 67, SiH	38, 67, SiH

Blue 1	Blue 2	Blue 3	Blue 4	Blue 5
				
38, 60, 67, SiH	38, 60, 67, SiH	38, 60, 67, SiH	38, 67, SiH	38, 67

Brown 1	Brown 2	Brown 3	Brown 4	Brown 5	Brown 6	Brown 7	Brown 8
							
38, 60, 67, SiH	38, 60, 67, SiH	38, 60, 67, SiH	38, 67, SiH	67, SiH	67, SiH	67, SiH	67, SiH

Grey 1	Grey 2	Grey 3	Grey 4	Grey 5	Grey 6
					
38, 60, 67, SiH	38, 60, 67, SiH	38, 60, 67, SiH	67, SiH	67, SiH	67, SiH

Violet 1	Violet 2	Violet 3
		
38, 60, 67	38, 60, 67	38, 67

Images are intended for representational purposes only

Iris and Pupil Sizes

Lens Diameter (mm)	Standard Iris Size (mm)
12.50	11.00
13.00	11.00
13.50	11.50
13.75	11.50
14.00	12.00
14.50	12.50
15.00	13.00
16.00	13.00
17.00	14.00
18.00	15.00

Available Iris Diameters (mm)

**10.00, 10.50, 11.00, 11.50, 12.00, 12.50, 13.00
13.50, 14.00, 14.50, 15.00, 15.50, 16.00, 16.50**

Clear / Open Pupil (mm)



2.00, 2.50, 3.00, 4.00, 4.50, 5.00, 6.00

Black / Closed Pupil (mm)



**2.50, 3.00, 3.50, 4.00, 4.50, 5.00, 5.50
6.00, 7.00, 8.00, 9.00**

Basic Tint



Basic tint colours may allow some of the natural iris colour through, this will be more noticeable with light tints on dark eyes

Opaque Tints



Opaque tints will darken the iris tint slightly and will not allow natural iris colours through

Recommended care system is Quattro multipurpose solution by Abatron, Hydrogen peroxide systems may be used. Other multipurpose solutions may cause tint to fade.

For any further technical advice please do not hesitate to call our Professional Services Team 01280 702002 Option 2